Supported Detectors	Conversion Screen	Communication
DRX-1 3543	Detached Gd2O2S:Tb	Wireless/Tethered
DRX-1C 3543	Deposited CsI (Tl)	Wireless/Tethered
DRX Plus 3543	Detached Gd2O2S:Tb	Wireless/Tethered
DRX Plus 3543C	Deposited CsI (Tl)	Wireless/Tethered
DRX Plus 4343	Detached Gd2O2S:Tb	Wireless/Tethered
DRX Plus 4343C	Deposited CsI (Tl)	Wireless/Tethered
DRX Plus 4343F	Detached Gd2O2S:Tb	Tethered only
DRX Plus 4343FC	Deposited CsI (Tl)	Tethered only
DRX Plus 2530C	Deposited CsI (Tl)	Wireless only
DRX-L	Detached Gd2O2S:Tb	Tethered only
LUX 35 (May need a MOD on	Deposited CsI (Tl)	Wireless/Tethered
the existing system)		
Focus HD 35	Deposited CsI (Tl)	Wireless only with charging
Focus HD 43	Deposited CsI (Tl)	Wireless only with charging

X-ray Tube - 1	
Type of x-ray tube	Varian RAD-60
Tube voltage range	40 to 150 kVp
Nominal focal spot size	0.6 and 1.2 mm focal spots
Target angle	12 degrees
Tube cooling rates	Anode cooling rate = 100,000 HU/minute
	Tube Housing cooling rate = 18,000 HU/minute
Heat storage capacity of anode	400,000 HU
Inherent filtration	0.7 mm Al at 75 kV

X-ray Tube - 2	
Type of x-ray tube	Varian RAD-92
Tube voltage range	40 to 150 kVp
Nominal focal spot size	0.6 and 1.2 mm focal spots
Target angle	12 degrees
Tube cooling rates	Anode cooling rate = 140,000 HU/minute
	Tube Housing cooling rate = 18,000 HU/minute
Heat storage capacity of anode	600,000 HU
Inherent filtration	0.7 mm Al at 75 kV

X-ray Tube – 3	
Type of x-ray tube	Canon XRR-4631G
Tube voltage range	40 to 150 kVp
Nominal focal spot size	0.6 and 1.2 mm focal spots
Target angle	12 degrees

Technical Specifications DRX-Evolution Plus SystemPublic Classification

AG6108 K

Uncontrolled unless otherwise indicated

Page: 1 of 12

Tube cooling rates	Maximum Anode Heat Dissipation=1690 HU/second
Heat storage capacity of anode	400,000 HU
Inherent filtration	1.1 mm Al at 75 kV

Collimator - 1	
Automatic collimator	Ralco R 221 DACS
	• Automatic collimation (PBL override capability)
	Manual knob control
	• Digital display of SID, longitudinal and transverse blade position,
	spectral filter selection (programmable to display inches or
	centimeters)
	6 pairs of lead faced shutter blades
Light output	Greater than 160 LUX at 100 cm (39.4 in.)
Lamp type	LED
X-ray proofing	< 40 mR/h with X-ray beam = 150 kVp/4 mA. EN 60601-1-3 par.
	29.204.3
Inherent filtration	2 mm Al at 75 kV
Additional spectral filtration	0 (none), 1 mm Al, 1 mm Al + 0.1 mm Cu, 1 mm Al + 0.2 mm Cu

Collimator - 2	
Manual collimator	Ralco R 221/A DHHS
	Manual knob control
	6 pairs of lead faced shutter blades
Light output	Greater than 160 LUX at 100 cm (39.4 in.)
Lamp type	LED
X-ray proofing	< 40 mR/h with X-ray beam = 150 kVp/4 mA. EN 60601-1-3 par.
	29.204.3
Inherent filtration	2 mm Al at 75 kV

Collimator - 3	
Automatic collimator	Ralco R 221 DACS (For Dual Energy)
	Automatic collimation (PBL override capability)
	Manual knob control
	Digital display of SID, longitudinal and transverse blade position,
	spectral filter selection (programmable to display inches or
	centimeters)
	6 pairs of lead faced shutter blades
Light output	Greater than 160 LUX at 100 cm (39.4 in.)
Lamp type	LED

Technical Specifications DRX-Evolution Plus SystemPublic Classification

AG6108 K

Uncontrolled unless otherwise indicated

X-ray proofing	< 40 mR/h with X-ray beam = 150 kVp/4 mA. EN 60601-1-3 par. 29.204.3
Inherent filtration	2.3 mm Al at 75 kV
Additional spectral filtration	0 (none), 1 mm Al+ 0.5 mm Ag , 1 mm Al + 0.1 mm Cu, 1 mm Al +
	0.2 mm Cu

Dose Area Product Meter (Optional Accessory)		
Туре	VacuTec VacuDAP 156 00 15	
Upper limit of response range	0.8 nC/uGy-m2	
Energy range	40 to 150 kV	
Energy response	-6% - 0% over the 50 to 150 kV range	
Chamber filter effect	0.2 mm Al (70 kV)	

X-ray Generator and PDU Asse	nbly	
High-frequency generator	Manufactured by Carestrea	m Health. Model CGN-80
PDU (power distribution unit)	Manufactured by Carestream Health	
Generator Output	• 65 kW high-frequency output with digital feedback control	
	circuitry	
	Optional 80 kW high-fre	quency output with digital feedback
	control circuitry	
	• 40 to 150 kVp	
	• 0.1 to 1000 mAs	
Range of generator output	<u>65kW</u>	<u>80kW</u>
	800 mA at 81 kVp	1000 mA at 80 kVp
	630 mA at 103 kVp	800 mA at 100 kVp
	500 mA at 130 kVp	630 mA at 126 kVp
	400 mA at 150 kVp	500 mA at 150 kVp
Power requirement	380/400/480VAC, 3-phase, 125kVA	
Generator and PDU size	Dimensions (H x W x D)	
	120.8 cm x 64.5 cm x 57.7	cm (47.6 in. x 25.4 in. x 22.7 in.)
Generator and PDU weight	245 kg (540 lb) Unpacked	

Technical Specifications DRX-Evolution Plus SystemPublic Classification

Uncontrolled unless otherwise indicated

Console	
Single console integrating	Software, image processing, and DICOM output by Carestream
exposure control, detector	Health
control, and image output	Computer manufactured by HP
	HP Engage Flex Pro (also called HP 5820) with i3 8100 3.6GHz (or
	better) or Z4 G4 with Intel XEON W-2225 4.1GHz (or better)
	Memory 16,32 GB
	Hard drive 512GB, 1TB
	10/100/1000 BaseT network interface
	DVD-RW drives
	Microsoft Windows 10
	Hard drive option: encrypted hard drive certified to meet U.S.
	National Institute of Standards Technology (NIST) Federal
	Information Processing Standard (FIPS) 140-2 Level 2
	Console PC specifications are subject to change by our workstation
	supplier. Workstation performance provided will be equal to or
	better than performance of the Standard or (optional) Enhanced PC
	below
Patient data entry	Touch-screen monitor; keyboard, bar code or DICOM work list
	through HIS/RIS (optional)
Selection of exposure	Full exposure parameters can be set manually or AEC (automatic
	exposure control) can be used. Anatomical programming can be set
	for a variety of examinations.
Image storage capacity	Up to 10,000 images can be saved for retransmission or
	reprocessing. Images can be protected; otherwise they are reclaimed
	on a FIFO (first in, first out) basis. Images can be saved to
	removable media on a DVD-RW drive.
DICOM 3.0 compliance	DICOM Work List: Comply (option);
	DICOM Store: Comply;
	DICOM Print: Comply;
	DICOM Modality Performed Procedure Step (option);
T 1.4 4 4	DICOM Store Commit: Comply
Image data output	12-bit log, 14.6 MB
Console Table (optional)	Dimensions (HxWxD): 93.3 cm x 90.0 cm x 60.0 cm $(36.7 \text{ in. x } 35.4 in.$
	1n. x 23.6 in.)
Other components	Interfaces with the detector and the x-ray exposure equipment while
	providing upstream and downstream connectivity to acquire and
	iransmit patient, exam, and image data in digital format.
	Detector Portable DRX-1, DRX- Plus
	• Preview image is available in <4 seconds (V5.5 SW or later);
	total cycle time (processing speed) is 20 seconds.

Technical Specifications DRX-Evolution Plus SystemPublic Classification

AG6108 K

Uncontrolled unless otherwise indicated

Console	
	 Detector Fixed Preview image is available in <3 seconds (V5.5 SW or later); total cycle time (processing speed) is 11 seconds.
	 Acquisition and distribution of data Operator input: log on/log off; patient demographics (add / edit patients); patient accession number, study selection; search by patient name, ID, or visit number; output setup selection (workstations, archives, printers); RIS interface to acquire patient demographic and exam data (optional); and auto transmission of patient/exam/image data.
	 Exposure and image capture Exposure factors for each exam view (programmed default factors with manual overrides; small/medium/large patient size selection; manual technique selection; tube warm-up capability; detector calibration). Automatic exposure control (AEC) (ion chamber selection; manual control, AEC on/off); Exam Tutor; detector array; detector rotation indicator (for wall stand only); exposure button (prep/exposure control button; audible/visible exposure indicator); and generator on/off.
	 Acquire and process digital images Optimize grayscale display and apply examination specific nonlinear edge enhancement with built-in perceptual tone-scale processing (PTS). Optional EVP (enhanced visualization image processing) software is available. Display preview image (apply image cropping; apply image multi-formatting; add image markers, flip and rotate image, enter technologist comments; accept/reject image).
	 Administrative System administrative (create/modify user logon and password; create/modify output printer configurations; manage image output queue; re-send image output; manage local database [view/delete patient image files]; protect selected patient/image files from reclamation; create/modify technologist comments; perform detector array calibration; generate test pattern images for output device QC). Remote diagnostic service capability via internet connection.

AG6108 K

Uncontrolled unless otherwise indicated

Supported image outputs
Class, DICOM Work List, and DICOM Print Service Class.

Monitor	
Туре	Multi-touch IPS panel with LED backlighting
Viewable Image Area	58.42 cm (23 in) wide screen diagonally measured
Display resolution	1920 x 1080, 60 Hz
Viewing angles	Typical 178° horizontal and vertical
Monitor Foot Print	With Stand 59.93 x 6.95 x 40.92 cm (23.6 x 2.7 x 16.1 in)
	Head only 59.93 x 4.05 x 50.7 cm (23.6 x 1.6 x 16.1 in)

Auto positioning and Motorizati	ion / Synchronization
Programmed Auto Positions	Automatic Positioning of the Overhead Tube Assembly, Table
(OPTIONAL FEATURE for	Bucky, and Wall Stand Bucky to any of 128 positions that are
Fully Automated)	programmed by the operator.
Motorization / Synchronization	Automatic Tracking of the Overhead Tube Assembly (X-ray beam
	axis) to follow the Table Bucky or Wall Stand Bucky position
	Automatic Centering and perpendicular alignment of the Overhead
	Tube Assembly (X-ray beam axis) to the Table Bucky or the Wall
	Stand Bucky position. (Fully Automated only)
	Automatic Centering of the Table Bucky to the Overhead Tube
	Assembly (X-ray beam axis) position.
	Automatic Centering of the Wall Stand Bucky Z axis (X axis with
	X-rail kit) to the Overhead Tube Assembly (X-ray beam axis)
	position. (Fully Automated and Hybrid Only)

Overhead Tube Assembly – Fully Automated	
Туре	5 axes of motion
	Manual movement in all 5 axes of motion,
	Motor Assist available in the x and y direction for a ramped
	acceleration as guided by the operator
	Motorized in all axis for auto tracking, auto centering, and auto
	positioning

Technical Specifications DRX-Evolution Plus SystemPublic Classification

Uncontrolled unless otherwise indicated

Range of movement	Tube rotation around the horizontal axis relative to zero (collimator pointing down) is +135 degrees and -180 degrees Tube rotation around the vertical axis relative to zero is +/- 180 degrees Range of telescope (vertical travel) is 161 cm (63.5 in.) Standard longitudinal rail length of 6 m (19.7 ft.) manual or motorized.
T II D O I	Transverse bridge length of 4 m (13.1ft.)
Tall Room Option	Telescope Extension, 15 cm
Controls	Single handle control and individual controls of longitudinal, transversal, vertical, and rotational directions; auto-centering and auto-tracking controls
Display	 Large 12.1 inch (diagonal) LCD touch-screen display, indicating: SID, Head Angle, Head Rotation Angle, mAs, kVp, Detector selection, Positioning status and Exam type. Allows toggling to Patient ID. Optional OTC Display Package to change techniques/views
Optional accessories	IR Remote

Wall Stand	
Range of Movement	• Up and Down: from floor to 180 cm
	Optional X-Rail: Up to 500 cm
	• Bucky Tilt: - 20 to 0 to 90 degrees (horizontal)
	• In horizontal position, Bucky center relative to column = 654 mm
	• In vertical position, Bucky surface relative to column = 605 mm
	Bucky rotation
	• Bucky Angulation L/R up to 45 deg in 15 deg increments
	(premium only)
Controls	• (Standard) – Vertical brake controls on both sides of wall stand
	• (Automatic) - Keypad (One Keypad is installed by default. A
	second keypad can be installed as an option) that can be
	configured/installed on both the premium wall stand, including
	motorized up/down and tilting as well as all collimator functions,
	auto positioning, auto centering, and auto tracking activation.
Motorization (Automatic)	Motorized motion for vertical, tilt, and optional rail drive
Automatic Exposure Control	• 5-cell ion chamber if equipped with DRX or DRX Plus detector
	• 3-cell ion chamber if equipped with non DRX detector
Bucky	The grid is stationary
	The grid is removable
	Grid present information provided to the console
	• Bucky housing assembly shaped for comfort in all positions; grids
	are removable in all convenient orientations (including under table)

Technical Specifications DRX-Evolution Plus SystemPublic Classification

AG6108 K

Uncontrolled unless otherwise indicated

Wall Stand	
Grids	 One grid provided; additional grids optional: General Purpose: 8:1, 40 lines per cm (103 lines per inch) focus140cm (55 in.). Interspace material: Aluminum Near: 12:1; 40 lines per cm (103 lines per inch) focus 110cm (43 in.). Interspace material: Aluminum Far & LLI (Different Grid Frames): 12:1; 40 lines per cm (103 lines per inch) Center focus 180cm (71 in.). Interspace material: Aluminum Far & LLI; 40 lines per cm (103 lines per inch) Center focus 180cm (71 in.) Interspace material: fibre Aluminum clad 3-slot floor standing grid holder (optional)
Grids (203L Option)	 Grid option: General Purpose: 12:1, 80 lines per cm (203 lines per inch) focus 140cm (55 in.) Near: 12:1; 80 lines per cm (203 lines per inch) focus 110cm (43 in.) Far & LLI (Different Grid Frames): 15:1; 80 lines per cm (203 lines per inch) focus 180cm (71 in.) Near-Pediatric: 10:1; 80 lines per cm (203 lines per inch) focus 110cm (43 in.) (recommended for pediatric imaging) Interspace material: aluminum Carbon fiber clad 3-slot floor standing grid holder (optional)
Accessories included	Lateral arm support
Optional accessories	BABIX holder

Fixed Elevating Table (VX3733-TABC)	
Туре	Elevating table, with four-way floating flat edge tabletop
Table movement	• Tabletop longitudinal travel: 122 cm (48 in.) range
	• Tabletop transverse travel: ± 13 cm (5.1 in.)
	• Table size 250 cm x 83 cm (98.4 in x 32.6 in)
	• Table height 52 cm (20.5 in.) minimum, 91 cm (35.8 in.) maximum.
	• Bucky longitudinal travel: 44 cm (17.3 in.)
	• Table maximum patient weight: 320 kg (705 lb)
	• Table maximum static load: 400 kg (881 lb)
Imageable area	• Width 63 cm (24.8 inch)
	• Length 250 cm (98.4 inch)
	Tabletop material: HPL
	• Radio-opacity is 0.7 mm aluminum equivalent at 100 kVp

Technical Specifications DRX-Evolution Plus SystemPublic Classification

AG6108 K

Uncontrolled unless otherwise indicated

Fixed Elevating Table (VX373)	3-TABC)
Motion controls	• Foot pedal operated elevation and tabletop float (double-step
	activation)
	Switch for automatic tracking of OTC longitudinally to Bucky
	• Switch for auto centering of OTC to Bucky
	Indicator for Grid Presence
	Option - Foot pedals located on rear of table
Automatic Exposure Control	3-cell ion chamber
Bucky	The grid is stationary
	The grid is removable
	Grid present information provided to the console
	• Detector assembly moves out of table for extremity views without
	tabletop (Grid remains in the table)
Grid	One grid provided:
	• 12:1, 40 lines per cm (103 lines per inch), focus 110 cm (43 in.),
	Interspace material: aluminum
	• 12:1, 40 lines per cm (103 lines per inch), focus 110 cm (43 in.),
	Interspace material: fiber
	Interspace material: aluminum
	Aluminum clad
	3-slot floor standing grid holder (optional)
Grid (203L Option)	Grid Option:
	• 12:1, 80 lines per cm (203 lines per inch), focus 110 cm (43 in.)
	Interspace material: aluminum
	Carbon fiber clad
	3-slot floor standing grid holder (optional)
Optional Accessories	Patient hand grips
	Compression strap
	Table pad
	Lateral cassette holder
	Tethered Remote

Fixed Elevating Table (VX3733-TAB)	
Туре	Elevating table with four-way floating tabletop
Table movement	• Tabletop longitudinal travel: 120 cm (47.2 in.) range
	• Tabletop transverse travel: ± 13 cm (5.1 in.)
	• Table size 240 cm x 83 cm (94.5 in x 32.5 in
	• Table height 53 cm (20.9 in.) minimum, 86 cm (33.8 in.)
	maximum.
	• Bucky longitudinal travel: 40 cm (15.7 in.)
	• Table maximum patient weight: 272 kg (600 lb)
Imageable area	• Width 63 cm (24.8 inch)

AG6108 K

Uncontrolled unless otherwise indicated

Fixed Elevating Table (VX3733	-TAB)
	• Length 202.5 cm (79.7 inch)
	Tabletop material: HPL
	• Radio-opacity is 0.7 mm aluminum equivalent at 100 kVp
Motion controls	• Foot pedal operated elevation and tabletop float (double-step
	activation)
	• Switch for automatic tracking of OTC longitudinally to Bucky
	Switch for auto centering of OTC to Bucky
	Indicator for Grid Presence
	Option - Foot pedals located on rear of table
Automatic Exposure Control	3-cell ion chamber
Bucky	• The grid is stationary
	• The grid is removable
	Grid present information provided to the console
	• Detector assembly moves out of table for extremity views without
	tabletop (Grid remains in the table)
Grid	One grid provided:
	• 12:1, 40 lines per cm (103 lines per inch), focus 110 cm (43 in.);
	Interspace material: aluminum
	• 12:1, 40 lines per cm (103 lines per inch), focus 110 cm (43 in.);
	Interspace material: fiber
	Aluminum clad
	3-slot floor standing grid holder (optional)
Grid (203L Option)	Grid Option:
	• 12:1, 80 lines per cm (203 lines per inch), focus 110 cm (43 in.)
	Interspace material: aluminum
	Carbon fiber clad
	3-slot floor standing grid holder (optional)
Optional Accessories	Patient hand grips
	Compression strap
	Table pad
	Lateral cassette holder
	Tethered Remote

Fixed Elevating Table (VX3733-TAB320)		
Туре	Elevating table, with four-way floating tabletop	
Table movement	• Tabletop longitudinal travel: 96 cm (37.8 in.) range	
	• Tabletop transverse travel: ± 13 cm (5.1 in.)	
	• Table size 216 cm x 83 cm (85.0 in x 32.5 in)	
	• Table height 53 cm (20.9 in.) minimum, 86 cm (33.8 in.) maximum.	
	• Bucky longitudinal travel: 40 cm (15.7 in.)	
	• Table maximum patient weight: 320 kg (705 lb)	

AG6108 K

Uncontrolled unless otherwise indicated

Fixed Elevating Table (VX3733-TAB320)		
Imageable area	• Width 63 cm (24.8 inch)	
	• Length 178.5 cm (70.3 inch)	
	Tabletop material: HPL	
	• Radio-opacity is 0.7 mm aluminum equivalent at 100 kVp	
Motion controls	• Foot pedal operated elevation and tabletop float (double-step	
	activation)	
	• Switch for automatic tracking of OTC longitudinally to Bucky	
	Switch for auto centering of OTC to Bucky	
	Indicator for Grid Presence	
	Option - Foot pedals located on rear of table	
Automatic Exposure Control	3-cell ion chamber	
Bucky	The grid is stationary	
	• The grid is removable	
	Grid present information provided to the console	
	• Detector assembly moves out of table for extremity views without	
	tabletop (Grid remains in the table)	
Grid	One grid provided:	
	• 12:1, 40 lines per cm (103 lines per inch), focus 110 cm (43 in.),	
	Interspace material: aluminum	
	• 12:1, 40 lines per cm (103 lines per inch), focus 110 cm (43 in.),	
	Interspace material: fiber	
	Interspace material: aluminum	
	Aluminum clad	
	3-slot floor standing grid holder (optional)	
Grid (203L Option)	Grid Option:	
	• 12:1, 80 lines per cm (203 lines per inch), focus 110 cm (43 in.)	
	Interspace material: aluminum	
	Carbon fiber clad	
	3-slot floor standing grid holder (optional)	
Optional Accessories	Patient hand grips	
	Compression strap	
	Table pad	
	Lateral cassette holder	
	Tethered Remote	

SMART (Optional Feature)	
Required Hardware	• Two camera modules (One camera attached to Collimator, other
(Please refer to full Smart	Lateral Position)
Feature specification AC2160	• Z4 G4 with Intel XEON W-2225 4.1GHz (or better) (Required PC,
for more details)	not included in Smart feature kit)
	• Full Motorized Overhead Tube Assembly with cable pre-installed

AG6108 K

Uncontrolled unless otherwise indicated

SMART (Optional Feature)	
	(Required, not included in Smart feature kit)
	• Z-motorized Wall Stand with Bucky Cover supporting Smart Feature
	(Required, not included in Smart feature kit)
Software Options	Smart Auto position
(Please refer to full Smart	Patient Position Monitoring
Feature specification AC2160	Smart Patient Positioning
for more details)	Smart Technique
	Smart Collimation
	Video Assistance
	Audio Assistance



"Rx only"

Technical Specifications DRX-Evolution Plus SystemPublic Classification

AG6108 K

Uncontrolled unless otherwise indicated

Page: 12 of 12