

X-DR Static Classic

Technical Specifications



Generator	32 kW / max. 450 mA
Voltage range	40 - 125 kV (optional throttling to 99 kV)
Tube	Monoblock, Focal spots 0.6 / 1.3 mm
Exposure time range	1 ms – 2.5 s
Programs	352 pre-programmed
Combi server	1 TB
Monitor	X-HW More Multitouch Monitor
Table top	160 cm x 75 cm (63" x 29.5"), floating Longitudinal ±35 cm (13.8") Transversal ± 11cm (4.3")
Mains connection	230 V - 50 / 60 Hz 1-phase (wall plug)
Detector	X-DR XL wired, 43 cm x 43 cm (17" x 17"), Csl
X-AQS PACS Vet	X-AQS PACS Vet X-AQS Module Archiving on CD/DVD X-AQS Module Data Compression contains for Client X-AQS Module Patient CD/DVD X-AQS Module Import Bitmap Images X-AQS Module Import DICOM Images X-AQS Module Import AQS Images X-AQS Module Image Export X-AQS Module Mail Export
X-AQS Client	1 workstation (expandable up to 4 workstations) X-AQS Viewer X-AQS Diagnosis License per Client, floating license for X-AQS Software